## **Answer Key**



## Lab Equipment Quiz

Wh	ich	is	larg	er?

2.	30	millig	rams	or <mark>30</mark>	grams	of suga	ır
3.	25	liters	or 25	millili	ters		

1. 5 kilograms or 5 grams of salt

Match the pieces of equipment to the descriptions:

4.	Used to measure volume of a liquidF	A. Thermometer
5.	Measures massD	B. Microscope
6.	Measures temperatureA	C. Meniscus
7.	Can be used to grow bacteriaG	D. Balance
8.	Used to see tiny thingsB	E. Tare
9.	The curve of found at the top of a liquidC	F. Graduated cylinder
		G. Petri dish

10. Describe the steps to measure the volume of an irregular object.

Water Displacement

- 1. Measure volume of water in a graduated cylinder
- 2. Add irregular object to graduated cylinder
- 3. Measure volume of water and irregular object in graduated cylinder
- 4. Subtract volume in step 3 from volume in step 1 to get volume of irregular object
- 11. How do you read the measurements on a triple-beam balance?

  Balance must be at 0. Add object. Adjust weighted drag slides on beams, starting with 100-gram beam, then 10-gram beam, and finally 1-gram beam, making sure the highest weight value is used from each beam without surpassing the object's weight until balance is once again centered at 0. Add weight measured on all three beams to obtain the object's total mass.